



“My husband and I feared he would get locked in and be unable to communicate... everyone with ALS is fearful of being locked in.




“It is important to have NeuroSwitch™ before there is a crisis... to accommodate his changing needs.



“Because of NeuroSwitch my husband can communicate with less effort and more speed than ever before...

...and there’s much less stress in our home as a result”*

Lisa, wife of USAF veteran Gary who has quadriplegia and loss of speech due to ALS



CONTROL  **SYSTEM**

NEURAL SYSTEM TECHNOLOGIES

BIONICS

TECHNOLOGIES

MIND OVER MATTER™

Control Bionics
333 North Broad Street
Fairborn, OH 45324
www.controlbionics.com
1(937)879-0734 (OFFICE)
1(202) 257-7090 (CELL)

NOTES:

* All testimonials are from NeuroSwitch clients or caregivers - details available on application

© River Tree Media 2009

LINSO-001F-20090915

Being
Locked in
Does Not
always
mean
being
Shut Out



Communicate Control Connect

NeuroSwitch users who can't use other Assistive Technology (AT) switches are:

- **controlling voice communicators**
- **surfing the web**
- **sending and receiving emails**
- **operating remotes for TVs & stereos**
- **setting environmental controls**
- **running regular computer programs**
- **able to run most other AT systems**

HOW NEUROSWITCH WORKS

NeuroSwitch™ uses signals from adhesive sensors on the skin - often in areas that are functionally paralyzed - to enable people with severe disabilities to control a switch. It uses proven and approved EMG monitoring with a patent-applied-for signal processing system.

Painless

Non-invasive

Intuitive

Master it in your first session

www.controlbionics.com

“My hands slipped off other AT switches...”

“...Once I'm connected to NeuroSwitch, I can forget about finding the switch... I become part of the system...”

“It's easier and smoother for me to use than any other switch.”*

Nick (Australia) ALS-Locked in Syndrome



“J'utilise avec bonheur neuroswitch tous les jours ... Le système fonctionne dès la mise sous tension et ne nécessite aucune retouche de la journée.” (“I am very happy using NeuroSwitch every day. The system operates as soon as you turn it on it isn't necessary to adjust it throughout the day.”)*

Emailed by Roger (France): physicist with ALS who designs electronic systems, surfs the web, and emails using NeuroSwitch to control his AT and design systems



SEE TV News story of Gary, and videos of others using NeuroSwitch at www.controlbionics.com

New in Spring 2010

NeuroSwitch G5 Liberator

DELIVERED FULLY INSTALLED ON A NEW

Apple MacBook Pro Laptop

ALL INCLUDED



“NeuroSwitch is now totally Mac OS X native, and wireless thanks to Bluetooth. It's fast, it's portable, it's punk rock cool.”*

Michael (USA): Apple Master, computer expert, author and blogger, uses NeuroSwitch to overcome Locked in Syndrome due to SMA. (Quote written using his NeuroSwitch).



VETERANS

NeuroSwitch has been approved for full cover by the VA for our first Gulf War veteran.

Contact us for information on assessments and applications.

www.controlbionics.com