




NEURONODE DUO ACCESS FULL POTENTIAL

The NeuroNode Duo is a speech generating device with powerful access methods that make the user more effective in communicating. The NeuroNode sensor allows full use of the communication software with the slightest movement.

Control Bionics offers a suite of dedicated speech generating devices and accessories to give people living with complex communication and physical needs their voice.

Our product range is more intuitive and flexible than ever.

CONTROLBIONICS.COM
+61 3 9897 3576



POWERFUL DEVICE

Established, trusted equipment to efficiently run the system.



EMG MODE

The slightest movement can activate the sensor.



SPATIAL MODE

A sensor that moves with the user to turn movement into selection.



COMMUNICATION SOFTWARE

Gives the user tools to express themselves and communicate independently.



FUNDING OPTIONS

Get the device your client needs, fully funded. We're a registered NDIS provider.



24/7 SUPPORT

You and your clients will always have a resource to help with the device.



POWERFUL ACCESS METHOD

The NeuroNode Duo features the NeuroNode, a sensor which offers the combined power of EMG and 3D spatial control. The NeuroNode is flexible and doesn't require specific body position or lighting conditions.

It allows users to choose muscle control or movement control to operate the NeuroNode Duo all in one small, wearable device.

With EMG (electromyography) and Spatial control options, the NeuroNode Duo gives the user communication possibilities they've never had before.

The wireless sensor moves with the user, allowing flexible placement for even the most minimal movement can trigger the sensor, giving the user control and independence.

STANDARD CONFIGURATION

The NeuroNode Duo is offered with a Windows Surface Pro tablet and a Control Bionics NeuroNode. This powerful combination is our standard device model and highly recommended.

OTHER CONFIGURATIONS

The NeuroNode Duo is available using an Apple iPad device with the same NeuroNode technology. Combining the flexible NeuroNode sensor, and an Apple iPad, this device lets the user operate their speech generating software using iOS.

WHAT'S IN THE BOX

- Windows Tablet
- Protective case
- Charger
- Mounting Frame
- VESA mounting bracket
- Communication software license
- NeuroNode device
- NeuroNode charger kit
- NeuroNode accessory kit
- NeuroNode Bands

Additional available options: Cloth, paediatric, facial, or snap electrodes



REGISTERED
NDIS
PROVIDER

