



POWERFUL ACCESS METHOD

The NeuroNode Duo features the NeuroNode, a sensor which offers the combined power of EMG and 3D spatial control. The NeuroNode is flexible and doesn't require specific body position or lighting conditions.

It allows users to choose muscle control or movement control to operate the NeuroNode Duo all in one small, wearable device.

With EMG (electromyography) and Spatial control options, the NeuroNode Duo gives the user communication possibilities they've never had before.

The wireless sensor moves with the user, allowing for flexible placement. Even the most minimal movement can trigger the sensor giving the user control and independence.

STANDARD CONFIGURATIONS

The NeuroNode Duo is offered with a Windows Surface Pro tablet and a Control Bionics NeuroNode. This powerful combination is our standard device model and highly recommended.

OTHER CONFIGURATIONS

The NeuroNode Duo is available using a Smartbox Grid Pad device with the same NeuroNode technology. The same great benefits found in the NeuroNode Duo, but now now available on a Smartbox Device.

The NeuroNode Duo is available using an Apple iPad device with the same NeuroNode technology. Combining the flexible NeuroNode sensor, and an Apple iPad, this device lets the user operate your speech generating software using iOS.

WHAT'S IN THE BOX

- Surface Pro Tablet
- Surface Pro Protective case
- Surface Pro Charger
- Surface Pro Mounting Frame
- Rehadapt VESA mounting bracket
- Grid 3 Full License
- NeuroNode device
- NeuroNode charger kit
- NeuroNode accessory kit
- NeuroNode Bands

Additional available options: Cloth, pediatric, facial, or snap electrodes

